

Cleanlant I-FIX System

I-FIX
system

DENTAL IMPLANT SOLUTION

C o n t e n t s

- A Multi-Center Retrospective Study of DENTIS Implant System3p
- Fixture System Selection4p
- Features of I-FIX5p
- Opening Sequence of Ampule8p
- I-FIX Post Type Fixture9p
- I-FIX Angled Type Fixture10p
- I-FIX O-Ring Type Fixture12p
- DENTIS Mini Implant Drilling Sequence13p
- DENTIS Mini Implant I-FIX KIT14p



I-FIX


A Multi-Center Retrospective Study of DENTIS Implant System

Thomas K. Lee, DDS
 Su-Kwan Kim, DDS, PhD
 Sang-Don Joo, DDS
 Sang-Chul Ko, DDS


Background: The design of DENTIS implant system, marketed since 2005 in Korea and abroad in 15 countries, is designed with the following characteristics: RBM surface treatments for time-proven osseointegration; tapered body with optimized thread designs for easy initial fixation at the time of placement surgery; 3 different abutment connection types for the same body design, allowing easier transition for the operator from existing systems in his/her armamentarium; and simplified prosthetic components. The purpose of this retrospective study was to evaluate clinical success rates for a new dental implant system called DENTIS in various private practice clinical settings.

Methods: 707 consecutive patients at 3 different clinical locations were treated with 1429 DENTIS implants. All 3 different abutment connection types of DENTIS Implant System (external hex connection type, octagonal conical-taper internal supra-gingival connection type, and submerged/bone level hex conical-taper internal connection type) were utilized in this study based on operator decision for each case. Implants were placed at various locations throughout the maxilla and mandible according to the treatment plan, including delayed and immediate placements after extractions. Various bone grafting procedures were done, including sinus augmentations, when clinically necessary. Patients were recalled and clinically examined at regular intervals along with radiographs to monitor clinical progression and prosthetic serviceability and stability.


Age Distribution



Recipients




Anatomic Locations




Results: Average time since implant placements were 26 months. Average time since delivery of prosthesis was 21 months. 27 implants out of 1429 implants had to be removed before delivery of definitive restorations for various clinical failure criteria, resulting in a failure rate of 1.9%. Cumulative survival rate was 98.1%. Average age of the patient population was 52 years old at the time of implant placement surgery, while youngest patient was 16 years old and oldest patient was 87 years old. 52.7% of the patient population was female, while 47.3% was male. While maxillary molar region had the highest risk of failures anatomically, diabetes and smoking were the highest medical condition risk factors. Prosthetic complication factors such as screw loosening, cemented crowns coming-off, and porcelain fractures affected 36 implants, resulting in 4.8% prosthetic complication rate for the 26 months of this study.


DENTIS Implant Success Rate




Failures




Failure Anatomic Sites



Failure Details / Factors



Prosthetic Complications




Conclusion: A new dental implant system called DENTIS is found to be performing well in various clinical practice settings in this retrospective multi-center study. Cumulative success rate of 98.1% for the DENTIS Implant System, as demonstrated in this study, compares favorably to most of the leading implant systems in the market now. Same patient population will continue be followed up in coming years for further evaluation of DENTIS implants.

Fixture System Selection

IFIX

Post Type


Collar: 2mm, 4mm
Length: 10mm, 13mm, 16mm



IFIX

Angled Type


Collar: 1.4mm, 3.4mm
Length: 10mm, 13mm, 16mm



IFIX

O-Ring Type

Collar: 1.5mm, 3.5mm
Length: 10mm, 13mm, 16mm



Features of I-FIX

1. Three Types of Fixture



Post Type

- Suitable for narrow space in the maxillary and mandibular anterior teeth
- One-body type design provides maximum strength for mastication




Post type Neck design

- Neck design helps to make fast healing time



Angle Type

- Specialized 2 pieces-design reinforce safe fixation
- Abutments are divided into Cemented and Angled for various cases and placement direction



Angle Type - Collar+Cuff

- Total gingival height is collar (1,4mm or 3,4mm) + cuff height(0,7)



O-Ring Type

- Suitable for the edentulous patients
- Creating fast and simple temporary prosthesis




O-Ring Type Collar

- Usable with regular size of retainer furthermore o-ring cap is available
- Collar heights are 1,5mm and 3,5mm

	Post Type	O-Ring Type	Angle Type	Angle Type	Cemented Abutment	Angle Abutment
Collar	2mm 4mm	1.5mm 3.5mm	1.4mm 3.4mm	→ → → Cuff	0.7mm	0.7mm

2. 2pieces-type with Maximum Strength


- 1,2 Hex Driver -



1. The fixation is strong and safe between screw and fixture with over 4 and half pitch screw connection
2. Morse taper connection between the implant and the abutment offers the friction locking and sealing.
3. Torque recommendation : 30-35Ncm


3. Another Function of Angled Type Mount

- Cover Screw & Healing Abutment




· Mount is used as healing abutment and cover screw after placement

4. Angle Type - Free Abutment



- Free, easy milling and prosthesis work as the application of DENTIS regular free abutment
- Groove in the abutment body - Strong fixation with cement.


5. For use where Narrow Interproximal Space



Impression Coping			
Height	12mm	18mm	23mm

Three types of impression coping available in various pick-up impression cases

18mm and 23mm impression copings provide convenient impression-taking in narrow interproximal space.




6. Easy Surgical Tool

One driver fits all 3 types of screw




Opening Sequence of Ampule

Post Type, O-Ring Type, Angle Type are the same way to open sequence of ampule.




1



2




3




4



5



6




7




8



9



A



B




● A tip for the installed implant fixture.

A. Remove the implant driver and fixture with a slight off-axis motion to disengage the connection before pulling away.

B. Do not set the drill torque to maximum torque of 80Ncm.

I-FIX Post Type Fixture

- Suitable for narrow space in the maxillary and mandibular anterior teeth
- One-body type design provides maximum strength for mastication
- Neck design helps to make fast healing time



Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.0	2.0mm	10.0mm	DMPF2010S	DMPF2010SN
	2.0mm	13.0mm	DMPF2013S	DMPF2013SN
	2.0mm	16.0mm	DMPF2016S	DMPF2016SN
	4.0mm	10.0mm	DMPF2010L	DMPF2010LN
	4.0mm	13.0mm	DMPF2013L	DMPF2013LN
	4.0mm	16.0mm	DMPF2016L	DMPF2016LN




Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.5	2.0mm	10.0mm	DMPF2510S	DMPF2510SN
	2.0mm	13.0mm	DMPF2513S	DMPF2513SN
	2.0mm	16.0mm	DMPF2516S	DMPF2516SN
	4.0mm	10.0mm	DMPF2510L	DMPF2510LN
	4.0mm	13.0mm	DMPF2513L	DMPF2513LN
	4.0mm	16.0mm	DMPF2516L	DMPF2516LN

Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø3.0	2.0mm	10.0mm	DMPF3010S	DMPF3010SN
	2.0mm	13.0mm	DMPF3013S	DMPF3013SN
	2.0mm	16.0mm	DMPF3016S	DMPF3016SN
	4.0mm	10.0mm	DMPF3010L	DMPF3010LN
	4.0mm	13.0mm	DMPF3013L	DMPF3013LN
	4.0mm	16.0mm	DMPF3016L	DMPF3016LN

Lab Analog	Code No.
	D MPL

I-FIX Angled Type Fixture

- Specialized 2 pieces-design reinforce safe fixation
- Abutments are divided into Cemented and Angled for various cases and placement direction
- Three types of impression coping available in various pick-up impression cases



Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.0	1.4mm	10.0mm	DMAF2010SM	DMAF2010SNM
	1.4mm	13.0mm	DMAF2013SM	DMAF2013SNM
	1.4mm	16.0mm	DMAF2016SM	DMAF2016SNM
	3.4mm	10.0mm	DMAF2010LM	DMAF2010LNM
	3.4mm	13.0mm	DMAF2013LM	DMAF2013LNM
	3.4mm	16.0mm	DMAF2016LM	DMAF2016LNM

Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.5	1.4mm	10.0mm	DMAF2510SM	DMAF2510SNM
	1.4mm	13.0mm	DMAF2513SM	DMAF2513SNM
	1.4mm	16.0mm	DMAF2516SM	DMAF2516SNM
	3.4mm	10.0mm	DMAF2510LM	DMAF2510LNM
	3.4mm	13.0mm	DMAF2513LM	DMAF2513LNM
	3.4mm	16.0mm	DMAF2516LM	DMAF2516LNM

Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø3.0	1.4mm	10.0mm	DMAF3010SM	DMAF3010SNM
	1.4mm	13.0mm	DMAF3013SM	DMAF3013SNM
	1.4mm	16.0mm	DMAF3016SM	DMAF3016SNM
	3.4mm	10.0mm	DMAF3010LM	DMAF3010LNM
	3.4mm	13.0mm	DMAF3013LM	DMAF3013LNM
	3.4mm	16.0mm	DMAF3016LM	DMAF3016LNM

Lab Analog	Code No.			
	DMAL			

Angled Abutment (15°)	Abutment Diameter	Cuff	Height	Code No.
	Ø4	0.7mm	9.1mm	DMAA15S


Angled Abutment Screw	Code No.			
	DMAS			

Cemented Abutment	Abutment Diameter	Cuff	Height	Code No.
	Ø4	0.7mm	8.3mm	DMCAHS
	Ø4	0.7mm	8.3mm	DMCANS

Free Abutment	Abutment Diameter	Cuff	Height	Code No.
 <p>0.7mm Hex Non-Hex</p>	<p>Ø4 Ø4</p>	<p>0.7mm 0.7mm</p>	<p>11.3mm 11.3mm</p>	<p>DMFAHS DMFANS</p>




Angled Mount Adapter	Detail	Code No.
	<p>0.9 Hex Driver</p>	<p>DMAMAS</p>

Gold UCLA	Abutment Diameter	Code No.
 <p>Hex Non-Hex</p>	<p>Ø4 Ø4</p>	<p>DMGUHS DMGUNS</p>

Impression coping	Type	Height(H)	Code No.
 <p>H 12mm H 18mm H 23mm</p>	<p>Pick-up</p>	<p>12mm 18mm 23mm</p>	<p>DMIC12S DMIC18S DMIC23S</p>

I-FIX O-Ring Type Fixture

- Suitable for the edentulous patients
- Creating fast and simple temporary prosthesis
- Usable with regular size of retainer furthermore o-ring cap is available



Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.0	1.5mm	10.0mm	DMOF2010S	DMOF2010SN
	1.5mm	13.0mm	DMOF2013S	DMOF2013SN
	1.5mm	16.0mm	DMOF2016S	DMOF2016SN
	3.5mm	10.0mm	DMOF2010L	DMOF2010LN
	3.5mm	13.0mm	DMOF2013L	DMOF2013LN
	3.5mm	16.0mm	DMOF2016L	DMOF2016LN

Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø2.5	1.5mm	10.0mm	DMOF2510S	DMOF2510SN
	1.5mm	13.0mm	DMOF2513S	DMOF2513SN
	1.5mm	16.0mm	DMOF2516S	DMOF2516SN
	3.5mm	10.0mm	DMOF2510L	DMOF2510LN
	3.5mm	13.0mm	DMOF2513L	DMOF2513LN
	3.5mm	16.0mm	DMOF2516L	DMOF2516LN


Implant Diameter	Collar	Length	Code	
			RBM	Machined
Ø3.0	1.5mm	10.0mm	DMOF3010S	DMOF3010SN
	1.5mm	13.0mm	DMOF3013S	DMOF3013SN
	1.5mm	16.0mm	DMOF3016S	DMOF3016SN
	3.5mm	10.0mm	DMOF3010L	DMOF3010LN
	3.5mm	13.0mm	DMOF3013L	DMOF3013LN
	3.5mm	16.0mm	DMOF3016L	DMOF3016LN

Lab Analog	Code No.
	DMOL

Cap	Code No.
	DOCS

O-ring	Code No.
	DOA0100 / DOA0400 / DOA0800

DENTIS Mini Implant Drilling Sequence



Ø2.0 Fixture	Point Drill → Ø1.3 Fixture Drill → Ø2.0 Fixture
Ø2.5 Fixture	Point Drill → Ø1.3 Fixture Drill → Ø1.8 Fixture Drill → Ø2.5 Fixture
Ø3.0 Fixture	Point Drill → Ø1.3 Straight Drill → Ø1.8 Straight Drill → Ø2.3 Straight Drill → Ø3.0 Fixture

DENTIS Mini Implant I-FIX KIT



DPD

Point Drill



DDE

Drill Extension



0,9 Hex
DRHDL09

Hex Driver




Ø1,3 Ø1,8 Ø2,3
DMD13 DMD18 DMD23

Straight Twist Drill




DRW

Ratchet Wrench



Code - DMK

Mini Machine Driver




Long Short
DMML DMMS

Mini Ratchet Driver



Long Short
DMRL DMRS

Depth Gauge



DMDP

Cleanlant

I FIX

www.dentisimplant.co.kr/eng/



Tel. +82-1899-2804 | Fax. +82-53-583-2806
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea



Tel. +1-323-677-4363~5 | Fax. +1-323-677-4366
6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +886-2-2808-5933 | Fax. +886-2-2624-2309
2F.,No.48,31,Minzu Rd., Tamsui Dist, New Taipei City 251,Taiwan(R.O.C.)



Tel. +021-5111-3828 | Fax. +021-5111-3828
上海市长宁区中山西路933号2205室