

Crestal Approach Surgical Procedure

OneQ Fixture Size Selection	Lindemann	Reamer Drill					Bone Condensing			Surgical Drill
	Ø2.0	Ø2.8	Ø3.1	Ø3.6	Ø4.1	Depth Gauge	Ovis Bone Allo/Xeno/BCP/HA	Bone Packer	OneQ Final Drill	
Ø3.9	●	●				●	○	●	Optional	
Ø4.2(4.1)	●	●	●			●	○	●		
Ø4.7(4.8)	●	●	●	●		●	○	●		
Ø5.2(6.0)	●	●	●	●	●	●	○	●		

● : Required ○: Optional

Lateral Approach Surgical Procedure

Technique Selection	Reamer Drill		Expanding	Sinus Elevation			Bone Graft
	LAD Ø8.0	LAR Ø8.0	Wall Expander	SAVE 01	SAVE 02	SAVE 03	Ovis Bone Allo/Xeno/BCP/HA
LAD 8.0 Technique	●		○	●	●	●	○
LAR 8.0 Technique		●	○	●	●	●	○

● : Required ○: Optional



SAVE
SINUS
Sinus Kits & Instruments



Tel. +82-1899-2804 | Fax. +82-53-583-2806
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea



Tel. +1-323-677-4363~5 | Fax. +1-323-677-4366
6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +021-5111-3828 | Fax. +021-5111-3828
上海市长宁区中山西路933号2205室

© 2018 DENTIS CO.,LTD. All Right Reserved.
Specifications are subject to change without notice.
Trademarks are the property of DENTIS CO., LTD.
Or their respective owners.

D-R-SAVE-V1-201811-EN



SAVE SINUS

SAVE SINUS KIT can be accessed from either crestal or lateral approach. It is designed to lift maxillary sinus membrane safer and faster. In addition, all components are specialized for DENTIS products and provide a successful implant solution when used with OneQ Fixture and Ovis Bone.

CRESTAL Approach



Lindemann Drill

- Drill Length : 13mm
- Diameter : Ø2.0
- Recommended RPM : 1,000RPM

|| CODE : DSSLD20



Crestal Stopper

- Prevents excessive progression of Crestal Reamer and instruments.
- Length : 3/4/5/6/7/8/9/10mm(8ea)



|| CODE : DSSCS3~10



Reamer

- Protect membrane perforation in case of Crestal Approach (Hatch Reamer)
- Crestal Stopper System (3~12mm) (Check bone thickness before Stopper attachment)
- Diameters : Ø2.8, Ø3.1, Ø3.6, Ø4.1
- Recommended RPM : 800~1,200RPM


|| CODE : DSSR28/31/36/41

CODE	DSSR28	DSSR31	DSSR36	DSSR41
Diameter(Ø)	Ø2.8	Ø3.1	Ø3.6	Ø4.1
SD Reamer				

Depth Gauge

- It is use to measure remaining bone thickness after SD Reamer.

|| CODE : DSSDG



Bone Packer

- It is to push bone material into drilled hole after SD Reamer.
- Crestal Stopper System (3~12mm)
- Caution** If remaining bone amount is little, refrain from using.
- Diameters : Ø3.5, Ø4.5
- After using SD-Reamer, It will be able to use for Sensor Gauge. (Can be used for membrane surface check)

|| CODE : DSSBP



Aqua Tip

- This tip is used to inject water when lift membrane. (Crestal Approach)

|| CODE: DSSAT

LATERAL Approach



Carbide Round Burr

- Burr for window open in Lateral Approach,
- Used for Low speed
- Diameter : Ø1.4, Ø1.8

|| CODE: DSSCB14L/18L



Diamond Round Burr

- Burr for window open in Lateral Approach.
- Used for Low speed
- Diameter : Ø2.0

|| CODE: DSSDB20L



LAD(Lateral Approach Drill)

- Curved blade relieve the impact of sinus membrane and generate bone lid.
- Lateral Stopper System (0.5~3.0mm) (Check bone wideness before Stopper attachment by image diagnosis)
- Diameters : Ø8.0((Recommend Wall Expander for enlarge window.)
- Recommended RPM : 800~1,200RPM

|| CODE: DSSLAD80



LAR(Lateral Approach Reamer)

- Blade design provides excellent cutting force and reduces the damage of the membrane by adopting a design in which the bonechip is filled in the blade during cutting.
- Lateral Stopper System (0.5~3.0mm) (Check image diagnosis bone width before stopper attachment.)
- Diameters : Ø8.0
- Recommended RPM : 800~1,200RPM

|| CODE: DSSAR80



Wall Expander

- Enlarge Lateral Window which performed by LAD or LAR.
- Crestal Stopper System(8/9/10/11/12) (Crestal Stopper Height=7mm=Lateral Wall width)
- Recommended RPM : 1,200RPM

|| CODE: DSSWE



Lateral Stopper

- Prevent excessive progress of Lateral Approach Drill/Reamer
- Length : 0.5/1.0/1.5/2.0/2.5/3.0mm(6ea)



|| CODE: DSSLS05~30

Sinus Curette

- To use for sinus membrane lift by its sequential order.
- Elevate Sinus Membrane accordance as SAVE 01/02/03 Instruments.
- To improve the activity in Maxillary sinus, tip size has been minimized.

|| CODE: DSSS01, 02, 03



Suction Tip

- Prevent membrane damage through proper suction force for Sinus surgery.

|| CODE: DSSST



Bone Well

- Used to store or mix bone/graft materials.

|| CODE: DSSBW