


# SQ-SL Fixture Specifications



A	Diameter	3.5	4.0	4.5	5.0	6.0	7.0	8.0
	Fixture							
B	Length	7 / 8 / 10 / 12 / 14				6 / 7 / 8 / 10 / 12		
C	Hex	1.7 Double Hex	2.5Hex					
D	Color							
	Guideline	Anterior	Anterior	Premolar	Molar	Molar	Molar	Molar


※ Length : 8~14mm - 0.5mm shorter than written length  
 ※ Nominal diameter does not correspond to the actual dimension

# SQ-SL Fixture Line-Up



Fixture Diameter	Length	Code No.
Ø3.5	7.0mm	DSSFN3507S
	8.0mm	DSSFN3508S
	10.0mm	DSSFN3510S
	12.0mm	DSSFN3512S
	14.0mm	DSSFN3514S

※ Set Code(S) : Fixture + Cover Screw[DSOFCS]



Fixture Diameter	Length	Code No.
Ø4.0	7.0mm	DSSFR4007S
	8.0mm	DSSFR4008S
	10.0mm	DSSFR4010S
	12.0mm	DSSFR4012S
	14.0mm	DSSFR4014S

※ Set Code(S) : Fixture + Cover Screw[DSCSM]



Fixture Diameter	Length	Code No.
Ø4.5	7.0mm	DSSFR4507S
	8.0mm	DSSFR4508S
	10.0mm	DSSFR4510S
	12.0mm	DSSFR4512S
	14.0mm	DSSFR4514S

※ Set Code(S) : Fixture + Cover Screw[DSCSM]